



 **InBody 270**






## Portable Analysis

Body composition analysis is key to helping your clients achieve their health or fitness goals, whether it's trying to lose fat or monitoring how much muscle is developing over time.

The InBody 270 goes beyond weight and measures **how much lean mass** you have **in each body segment**.

With these values, you can monitor how your body is adjusting to dietary changes, fitness routines, and lifestyle modifications. Its portable and foldable design is specifically designed for professionals who need to analyze body composition on-the-go.

### Features

-  **NO ESTIMATIONS**  
Only impedance is used to calculate your results; no statistical data needed
-  **15 SECONDS**  
Take a quick and easy body composition test
-  **HISTORY**  
Track progress with the body composition history chart on the results sheet
-  **LEAN MASS**  
See lean mass values for each body segment in pounds
-  **PORTABLE**  
Transport the unit in a carrying bag for analysis on-the-go



# Sample InBody 270 Results Sheet

# InBody

[InBody270]


SEE WHAT YOU'RE MADE OF

ID	Height	Age	Gender	Test Date / Time
Jane Doe	5ft.01.8in.	51	Female	05.04.2015 09 : 46


## Body Composition Analysis

Total amount of water in body	<b>Total Body Water</b> (lbs)	60.0
For building muscles and strengthening bones	<b>Dry Lean Mass</b> (lbs)	21.6
For storing excess energy	<b>Body Fat Mass</b> (lbs)	48.7
Sum of the above	<b>Weight</b> (lbs)	130.3


## Muscle-Fat Analysis




## Obesity Analysis



## Segmental Lean Analysis



## Body Composition History



## Body Fat-Lean Body Mass Control

**Body Fat Mass** - 22.0 lbs  
**Lean Body Mass** + 8.4 lbs  
 (+) means to gain fat/lean (-) means to lose fat/lean

**Lean Body Mass** 81.6 lbs

**Basal Metabolic Rate** 1168 kcal

## Results Interpretation

### Body Composition Analysis

The body weight is the sum of Body Fat Mass and Lean Body Mass, which is composed of Dry Lean Mass and Total Body Water.

### Muscle-Fat Analysis

Compare the bar lengths of Skeletal Muscle Mass and Body Fat Mass. The longer the Skeletal Muscle Mass bar is compared to the Body Fat Mass bar, the stronger the body is.

### Obesity Analysis

BMI is an index used to determine obesity by using height and weight. PBF is the percentage of body fat compared to body weight.

### Segmental Lean Analysis

Evaluates whether the amount of muscle is adequately distributed throughout the body. Compares muscle mass to the ideal.

### Body Composition History

Track the history of the body compositional change. Take the InBody Test periodically to monitor your progress.

### Body Fat-Lean Body Mass Control

Based on current body composition, the recommended change in Lean Body Mass and Body Fat Mass for a good balanced ratio. The '+' means to gain and the '-' means to lose.

### Basal Metabolic Rate

Basal Metabolic Rate is the minimum number of calories needed to sustain life at a resting state. BMR is directly correlated with Lean Body Mass.

## Results Interpretation QR Code

Scan the QR Code to see results interpretation in more detail.



## Impedance

Z(Ω)20 kHz	RA	LA	TR	RL	LL
100 kHz	345.0	358.5	23.4	286.6	296.0
	322.0	335.5	21.2	273.2	282.6



Copyright ©1996- by InBody Co., Ltd. All rights reserved. BR-USA-F9-A-141128

### Frequencies

20, 100 kHz

### Test Duration

15 seconds

### Dimensions

14.0 x 31.3 x 39.1 (L x W x H) : in

### Equipment Weight

30.9 lbs

### Database

100,000 results (if member ID is utilized)

### Warranty

1 Year Manufacturer's Warranty

### Weight Range

22-551 lbs

### Age Range

3-99 years

### Height Range

3 ft 1.4 in-7 ft 2.6 in

### Compatible Printers

Laser/Inkjet PCL 3 or above, SPL

### Measurements

10 impedance measurements  
 2 frequencies at each of the 5 segments  
 (Right Arm, Left Arm, Trunk, Right Leg, Left Leg)

### Additional Features

Lookin'Body 120 and Lookin'Body Web Compatible, Touch Screen, Voice Guidance System, Wi-Fi/Bluetooth Connectivity, Security Access Code, Foldable Design

### Accessories

InBody Results Sheets, InBody Tissues, USB Thumb Drive, Carrying Case, Thermal Printer

### Outputs

#### Standard outputs

Weight, Total Body Water, Dry Lean Mass, Lean Body Mass, Body Fat Mass, Skeletal Muscle Mass, Body Mass Index, Percent Body Fat, Segmental Lean Analysis, Body Composition History, Body Fat-Lean Body Mass Control, Basal Metabolic Rate, Segmental Impedance at each Frequency